

## Assessment of twelve echovirus virus-neutralisation assays in Europe: recommendations for harmonisation of non-polio enterovirus sero-surveillance studies

Karen Couderé, Kimberley S.M. Benschop, Gerrit Koen, Hetty van Eijk, Heli Harvala, Jean-Luc Bailly, Maël Bessaud, Everlyn Kamau, Isabelle Simon, Marie-Line Joffret, et al.

## ► To cite this version:

Karen Couderé, Kimberley S.M. Benschop, Gerrit Koen, Hetty van Eijk, Heli Harvala, et al.. Assessment of twelve echovirus virus-neutralisation assays in Europe: recommendations for harmonisation of non-polio enterovirus sero-surveillance studies. Journal of General Virology, 2024, 105 (9), 10.1099/jgv.0.001986. pasteur-04711395

## HAL Id: pasteur-04711395 https://pasteur.hal.science/pasteur-04711395v1

Submitted on 27 Sep 2024

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

# Assessment of twelve echovirus virus-neutralisation assays in Europe: recommendations for

### 2 harmonisation of non-polio enterovirus sero-surveillance studies

3	Karen Couderé <sup>1,2,*</sup> , Kimberley S.M. Benschop <sup>3</sup> , Gerrit Koen <sup>2</sup> , Hetty van Eijk <sup>2</sup> , Heli Harvala <sup>4,5</sup> , Jean-Luc
4	Bailly <sup>6</sup> , Maël Bessaud <sup>7</sup> , Everlyn Kamau <sup>8</sup> , Isabelle Simon <sup>6</sup> , Marie-Line Joffret <sup>7</sup> , Lubomira Nikolaeva-Glomb <sup>9</sup> ,
5	Irina Georgieva $^9$ , Asya Stoyanova $^9$ , Sabine Diedrich $^{10}$ , Sindy Böttcher $^{10}$ , Maria Cabrerizo $^{11}$ , Irena
6	Tabain <sup>12</sup> , Željka Hruškar <sup>12</sup> , Vladimir Stevanović <sup>13</sup> , Petri Susi <sup>14</sup> , Eero Hietanen <sup>14,15</sup> , Paula Palminha <sup>16</sup> ,
7	Petra Rainetová $^{17}$ , Anda Baicus $^{18}$ , Maartje Kristensen $^1$ , Dasja Pajkrt $^{19}$ , Jean-Luc Murk $^1$ , Katja C.
8	Wolthers <sup>2</sup> and European non-polio enterovirus network <sup>20</sup>
9	
10	
11	$^1$ Microvida, Laboratory of Medical Microbiology and Immunology, Elisabeth-TweeSteden Hospital, Tilburg,
12	Netherlands; <sup>2</sup> Department of Medical Microbiology, OrganoVIR Labs, Amsterdam UMC, AMC,
13	Amsterdam, Netherlands; <sup>3</sup> National Institute for Public Health and the Environment, RIVM, Bilthoven,
14	Netherlands; <sup>4</sup> Microbiology Services, NHS Blood and Transplant, Colindale, UK; <sup>5</sup> Radcliffe Department
15	of Medicine, University of Oxford, Oxford, UK; <sup>6</sup> Université Clermont Auvergne, CNRS, LMGE, UFR de
16	Médecine et des Professions Paramédicales, Clermont-Ferrand, France; <sup>7</sup> Institut Pasteur, Université de
17	Paris Cité, CNRS UMR 3569, Virus Sensing and Signaling Unit, Laboratoire associé au Centre national de
18	référence Entérovirus/Paréchovirus, Paris, France; <sup>8</sup> Nuffield Department of Medicine, University of Oxford,
19	Oxford, UK; <sup>9</sup> Department of Virology, National Center for Infectious and Parasitic Diseases, Sofia,
20	Bulgaria; <sup>10</sup> Robert Koch Institute, Berlin, Germany; <sup>11</sup> Enterovirus and Viral Gastroenteritis Lab, National
21	Centre of Microbiology, Instituto de Salud Carlos III, Consortium of Epidemiology and Public Health

22	(CIBERESP), Translational Research Network in Pediatric Infectious Diseases (IdiPaz), Madrid, Spain;
23	<sup>12</sup> Croatian Institute of Public Health, Zagreb, Croatia; <sup>13</sup> Faculty of Veterinary Medicine, Zagreb, Croatia;
24	$^{14}$ Institute of Biomedicine, Faculty of Medicine, University of Turku, Turku, Finland; $^{15}$ School of
25	Molecular and Cellular Biology, Faculty of Biological Sciences, University of Leeds, Leeds, UK; <sup>16</sup> Nacional de
26	Referência das Doenças Evitáveis pela Vacinação, National institute of Health, Lisbon, Portugal; <sup>17</sup> National
27	Institute Public Health NRL Enteroviruses, Praha, Czech Republic; <sup>18</sup> Cantacuzino National Institute for
28	Medical-Military Research and Development, Bucharest, Romania; <sup>19</sup> Emma Children's Hospital,
29	Department of Pediatric Infectious Diseases, OrganoVIR Labs, Amsterdam UMC, AMC, University of
30	Amsterdam, Amsterdam, Netherlands; <sup>20</sup> https://escv.eu/european-non-polio-enterovirus-network-
31	enpen/.
32	
33	*Correspondence: Karen Couderé, k.coudere@etz.nl
34	
35	
36	Keywords: Virus neutralisation assay; enterovirus; echovirus; surveillance; standardisation.
37	
38	Abbreviations: A549, alveolar carcinoma epithelial cell line; AFM, acute flaccid myelitis; Caco-2,
39	colorectal adenocarcinoma cell line; CPE, cytopathogenic effect; CVA, coxsackievirus A; CVB, coxsackievirus
40	B; DMSO, dimethylsulfoxide; E, echovirus; ELISA, enzyme-linked immunosorbent assay; EMEM, eagle's
41	minimum essential medium; ENPEN, European non-polio enterovirus network; EV, enterovirus; GCV,
42	geometric coefficient of variation; HT29, human B cell lymphoma cell line; IVIG, intravenous
43	immunoglobulins; LLC-MK, rhesus monkey kidney epithelial cell line 2; nAb, neutralising antibodies;

44 NPEV, non-polio enterovirus; RD, rhabdomyosarcoma; SOP, standard operating procedure; TCID50, 50% tissue culture infectious dose; VNA, virus neutralisation assay.

Abstract 46

45

Non-polio enteroviruses (NPEV) cause significant disease worldwide. Population-based sero-surveillance, 47 48 by measuring anti-bodies against specific NPEV types, provides additional information on past 49 circulation and the prediction for future upsurges. Virus neutralisation assays (VNA), the current method of choice for measuring NPEV type specific antibodies, are not entirely standardised. Via the European 50 Non-Polio Enterovirus Network, we organised a VNA quality assessment in which twelve labo- ratories 51 participated. We provided five echovirus (E) types (E1, E18, E30 G2, E30 G6 and E6) and intravenous 52 53 immunoglobulins (IVIG) as a sample for the NPEV VNA quality assessment. Differences in VNA protocols 54 and neutralising Ab (nAb) titres were found between the participating laboratories with geometric coefficients of variation ranging from 10.3–62.9%. Mixed-effects regression analysis indicated a small 55 but significant effect of type of cell line used. Harmonisation of cell line passage number, however, did 56 57 not improve variation between laboratories. Calibration by making use of a reference sample, reduced variation between laboratories but differences in nAb titres remained higher than two log2 dilution steps. 58 59 In conclusion, sero-surveillance data from different laboratories should be compared with caution and standardised protocols are needed. 60

- 61
- 62
- 63
- 64
- 65

#### 67

#### 68 INTRODUCTION

Non-polio enteroviruses (NPEV) are a significant source of disease worldwide [1, 2]. More than 110 69 70 NPEV types exist, including coxsackieviruses A and B (CVA and CVB), echoviruses (E) and the 71 numbered enteroviruses (EV). They are the leading cause of viral meningitis in children [3, 4]. Types like EV-D68, EV-A71, E11 and E30 have caused important upsurges of acute flaccid myelitis (AFM), 72 73 meningoencephalitis and myocarditis respectively [5–9]. The European Non-Polio Enterovirus Network (ENPEN) brings together different stakeholders across Europe to improve surveillance and research of 74 NPEV [10, 11]. In general, molecular surveillance is performed by genotyping EV positive samples 75 76 obtained from symptomatic cases. As the circulation of NPEV is mainly driven by population 77 immunity, seroprevalence data can help to predict upsurges and the susceptibility of a population to 78 a new viral strain [12]. As ELISA methods suffer from problems like cross-reactivity that do not 79 necessarily reflect immunity, virus neutralisation assays (VNA) are the method of choice to define 80 population immunity. Sero-surveillance studies using this method have been published by different 81 research groups [13–18]. However, measuring type-specific antibodies to the more than 110 NPEV 82 types is complex and cross-neutralisation has been reported [19]. Furthermore, VNA procedures are not entirely standardised, which could influence assay neutralising Ab (nAb) titres and hamper 83 comparison of data. With this study, we assess capacity for standardised sero-surveillance studies 84 85 within the ENPEN network.

#### 88 METHODS

89 Survey VNA

90 An invitation to participate in the study was sent out to different ENPEN laboratories in April 2021.

- 91 Twelve laboratories from Bulgaria, the Czech Republic, Finland, France (n=2), Germany, Croatia, the
- 92 Netherlands, Portugal, Spain, Romania, and the United Kingdom (England) agreed to participate in
- 93 the first quality assessment reported here. Participating laboratories were asked to provide their VNA
- 94 standard operating procedure (SOP) (summarised in Table 1).
- 95

#### 96 First quality assessment

97 Five echovirus types E1 (Virus [V]1), E18 (V2), E30 G2 (genogroup II) (V3), E30 G6 (genogroup III) (V4) 98 and E6 (V5) (Table S1, available in the online version of this article) and four individual aliquots of IVIG (Nanogam, lot number 17D29H462A, 2017, Prothya Biosolutions, the Netherlands) dilutions 1:1 99 (Serum [S]1), 1:4 (S2), 1:16 (S3), 1:64 (S4) were sent at room temperature to the participating 100 laboratories using the V and S codes. In case the estimated transport time took >48 h (laboratories 3, 101 4 and 12), transportation was done on dry ice. An instruction form was provided including the 102 103 guideline to store virus stocks at -80°C and samples at -20°C until performance of the assay. 104 Participating laboratories were asked to determine the 50% tissue culture infectious dose (TCID<sub>50</sub>) ml<sup>-1</sup> of V1–V5 and nAb titres according to the laboratories' own VNA SOP against V1–V5 in S1 and 105 against V5 in S2–4. The nAb results were collected on a standard form (added as File S1). Inter-106 laboratory variation of nAb titres was analysed by calculating geometric coefficients of variation (CV) 107

of the log2 transformed nAb titres against V1–V5 in S1. To assess intra-laboratory variation and to
 evaluate the possibility to use IVIG as a reference sample the linearity of the provided nAb titres
 against V5 in the IVIG dilution series was analysed for each laboratory by calculating Spearman's rank
 correlation coefficients.

112

#### 113 Second quality assessment: harmonisation and calibration

The second quality assessment aimed to reduce inter-laboratory variability by introducing 114 reference samples, to be used as calibrators, and by standardising the cell line (RD99) and day of 115 read-out. Ten of the twelve initial laboratories from Bulgaria, Finland, France (n=2), Germany, the 116 117 Netherlands, Portugal, Spain, Romania, and the United Kingdom (England) participated in the second assessment. Two echovirus types (E30 G2 [V3] and E30 G6 [V4]) and three IVIG dilutions 1:1 (S1), 1: 118 119 4 (S2), 1:16 (S3) were sent to the participating laboratories. A reference serum (R) dilution series 1:8 120 (R1), 1:16 (R2), 1:32 (R3), and 1:64 (R4) and a RD cell line (RD99 cell line passage 235) frozen in 121 dimethylsulfoxide (DMSO) 10% and foetal bovine serum 10% was also provided. All materials were sent on dry ice. Three of the participating laboratories (2, 10 and 11) experienced problems with 122 culturing or contamination of the frozen cell line stock. For these laboratories a flask with a 123 124 monolayer in Eagle's Minimum Essential Medium (EMEM) and FBS 8% was sent at room 125 temperature. An instruction form was provided including the guideline to store virus stocks and cell lines at  $-80^{\circ}$ C and samples at  $-20^{\circ}$ C until performance of the assay. We asked to determine nAb 126 titres in S1–S3 against V3 and V4 according to the laboratories' VNA SOP with the provided RD99 and 127 laboratories' proper cell lines and to perform an additional CPE read-out at day 7. We asked to 128 129 determine nAb titres in R1–R4 according to the laboratories' VNA SOP with the provided RD99 cell

130	line, and to perform an additional CPE read-out at day 7. As we did not know the real nAb titre of
131	the reference sample R, we took the median of measured nAb titres in R1 of all participating
132	laboratories and we calculated nAb titres in R2–R4 by considering their dilution factors. For each
133	laboratory, the measured log2 transformed nAb titres in S1–S3 were calibrated according to the
134	regression formula calculated based on R1–R4 (Fig. S1). By calculating the geometric CV values, the
135	dispersion was compared between measured (non-calibrated) and calibrated nAb titres in S1–S3.
136	
137	Mixed-effects regression analysis
138	To evaluate the individual effects of different VNA variables on the nAb titres, we performed a mixed-
139	effects regression analysis. Log2-transformed nAb titres were predicted by cell line, initial TCID <sub>50</sub> , virus
140	inactivation, virus/sample co-incubation (hours) and day of read-out as fixed effect independent
141	variables. To avoid overfitting of the model different cell line passages were not considered in the
142	model. Cell lines were divided into two categories: standardised cell line (RD99) or other cell lines.
143	Virus and laboratory were added as random effects. The significance of each independent fixed-effect
144	variable was calculated by comparing the likelihood of the model with and without the respective
145	variable via an ANOVA likelihood ratio test.
146	
147	Figures and statistics

All figures and statistics were conducted in R version 4.0.3 in R studio version 2023.03.1+446 with
packages ggplot2, dplyr, PKNCA and Ime4 [20].

150

151 **RESULTS** 

#### 152 VNA characteristics

Similarities and differences in VNA characteristics are presented in Table 1. All laboratories 153 except laboratory seven co-incubated 100 TCID<sub>50</sub> of the virus with sample dilutions, incubation 154 lengths ranging from 1 to 3 h. More than half of the participants performed serum sample 155 inactivation (7/12). There were differences in type of cell lines used, whether they were used in 156 157 suspension or monolayer and number of cells added. Nine of the twelve participating laboratories used RD as a cell line for all included echovirus types. Seven of the twelve participating laboratories 158 used a monolayer of cells. The number of cells, when counted, ranged from 14400 to 25000 cells per 159 well. Assay incubation length ranged from 3 to 7 days. All participating laboratories used CPE or 160 Crystal Violet (as an indication of CPE absence) as a read-out. Two main methods of nAb titre 161 calculation were used: the last dilution with 100% neutralization (both wells with less than 50% CPE) 162 (n=4) and the Reed Muench method (n=7) [21]. 163

164

#### 165 First quality assessment

166 The first quality assessment aimed to evaluate the inter-laboratory variability of VNA titres.

167 Differences of approximately three log10 TCID<sub>50</sub> ml<sup>-1</sup> were observed when culturing the echovirus

- 168 types in the laboratories' conditions (Table S2). A wide range of nAb titres was observed between
- 169 laboratories, especially against E1 and E18 (Fig. 1a and Table S3). Neutralising Ab titres ranged from 4 to
- 170 >1024 for E1 and E18, from 6 to >1024 for E30 G2, from 8 to 1024 for E30 G6 and from 64 to 1024 for E6.
- 171 The geometric CV for E1, E18, E30 G2, E30 G6 and E6 were 51, 62.9, 40, 31.7 and 17.2% respectively. Five
- 172 of the twelve participating laboratories reported back titration results, of which the majority were

between  $10^{1.5}$  and  $10^{2.5}$  TCID<sub>50</sub> per well (Table S3). Spearman's rank correlation coefficients of nAb titres in the IVIG dilution series, reported by nine of the twelve participating laboratories, ranged from -0.95 to -1 (Fig. 1b)

176

177 Second quality assessment

Our second quality assessment aimed to reduce inter-laboratory variability by introducing reference 178 samples to be used as calibrators and by standardizing the cell line (RD99) and day of read-out. 179 Neutralising Ab titres against E30 G2 ranged from 256 to 1024 and 256 to 2048 with the respective 180 181 laboratories' proper and RD99 cell lines, corresponding to a geometric CV of 8.2 and 10.5% respectively 182 (Fig. S2). Neutralising Ab titres against E30 G6 ranged from 128 to 640 and 128 to 1024 with the 183 respective laboratories' proper and RD99 cell lines, corresponding to a geometric CV of 13.1 and 184 13.9%, respectively. Although a large spread of initial TCID<sub>50</sub> results between different laboratories 185 was observed, there were no major intra-laboratory differences in initial TCID<sub>50</sub> results between the provided RD99 and the laboratories' cell lines (mostly RD cell lines) (Table S3). An exception was 186 laboratory nine, where the TCID<sub>50</sub> ml<sup>-1</sup> of E30 G2 was tenfold higher with the Caco-2 cell line than 187 188 with the RD99 cell line.

189 Back titration results ranged from 4 to 3162 TCID<sub>50</sub> per well (Table S4). A very weak correlation was

found between nAb titres and  $TCID_{50}$  back titration results (r=-0.18) (ig. S3). There were no intra-

191 laboratory differences in nAb titres between different days of read-out, ranging from 3 to 7 days.

192 Spearman's rank correlation coefficients in R1–R4 ranged from 0.8 to 1. Calibration, according to the

regression formula calculated based on R1–R4, reduced the dispersion of nAb titres in S1–S3 with in

194 general lower geometric CVs between laboratories for calibrated nAb titres compared to measured

195 nAb titres. Several outliers remained (Fig. 2 and Table 2).

196

- 197 Mixed-effects regression analysis
- 198 Since no clear cause of the nAb titre variation could be identified after the first and second quality
- assessment, we analysed the independent effects of different VNA variables (cell line, initial TCID<sub>50</sub>,
- 200 virus inactivation, virus/sample co-incubation [hours]
- and day of read-out), controlling for the random effect of virus (E1, E18, E30 G2, E30 G6 or E6) and
- 202 participating laboratory. Cell line (RD99 or other) independently affected nAb titre (P=0.02651),
- lowering the nAb titre by 0.804190 log2 or 1.75 when using RD99. Although statistically significant,
- this is less than what is considered normal variation (one log2 dilution step). No other significant
- 205 effects were observed.

206

#### 207 DISCUSSION

208 In this study, we tried to investigate the comparability of sero-surveillance data by organizing a VNA guality assessment and comparing different SOPs of twelve laboratories of the ENPEN network 209 across Europe. A wide variation of nAb titres was found, with differences ranging from 16- to 256-210 fold. For two of the five viruses tested (E1 and E18), nAb titres ranged from negative (below cut-off) 211 to >1024 (above cut-off). Along with differences in nAb titres, there were differences in VNA 212 213 characteristics like sample inactivation, incubation time, use of cell line (type, suspension/monolayer 214 and number of cells), day of read-out and nAb titre calculation. This means that results from serosurveillance studies between different laboratories cannot be compared for these echoviruses, and 215

216 most probably for other NPEV neither. Inter-laboratory differences between VNA nAb titres can be improved using a reference or an international standard serum sample [22]. In our study, inter-217 laboratory variation improved after using a type-specific animal serum sample as a calibrator but 218 remained considerable. It is possible that non-systematic errors, such as calculation or pipetting 219 errors and differences in read-out interpretation caused some of the outliers. VNA are laborious 220 assays susceptible to inter-individual variation. Further automatisation and standardisation for 221 example of read-out techniques (measurement of viral antigens or cell death using ELISA) could 222 potentially reduce inter-assay variation [23]. With more than 110 enterovirus types, providing a 223 reference sample for every individual type would be a major challenge. Therefore, we evaluated if 224 IVIG could serve as a universal NPEV reference sample. We found good linearity against all analysed 225 226 echovirus types, making IVIG a possible and more practical alternative to type-specific reference 227 serum samples. Importantly, IVIG of one identical batch should be used. Associations between nAb 228 titres and different VNA variables were explored. High differences in echovirus TCID<sub>50</sub> results between different laboratories were found, and preliminary analysis suggested a possible association 229 between nAb titres and cell line passage, day of read- out and back titration results. We therefore 230 standardised the cell line and day of read-out. Mixed-effects regression analysis indicated a small but 231 232 significant effect of type of cell line used. Our hypothesis was that differences in TCID<sub>50</sub> results for the 233 same quantity of virus are caused by differences in cell line susceptibility. This results in differences of virus added to the VNA and thus resulting in differences in nAb titres [24]. In the second assessment 234 however, TCID<sub>50</sub> results did not differ between the laboratories' RD cell lines and the provided RD99 235 cell line. Moreover, the dispersion of nAb titres was not reduced by harmonising the cell line. We did 236 see a substantial difference between the TCID<sub>50</sub> results on Caco-2 and RD99. Accordingly, the cell line 237 238 type probably plays a more important role than the number of cell line passage. As most laboratories

239 used RD cell lines it could be that we underestimated the effect of harmonising the cell line. Furthermore, other unknown factors in the treatment of these cell lines could play a role. All 240 possible factors causing variation need to be addressed. While awaiting standardisation via 241 international agreement on assay procedures and references samples (like for poliovirus and EV-A71) 242 [25, 26] laboratories could, when performing sero-surveillance studies, send a subset of the samples 243 and the respective virus strain to another laboratory for confirmation. Or if several laboratories are 244 involved in a sero-surveillance study, the same EV strain and a reference sample should be agreed 245 246 upon.

There are several limitations to our study. One is that only a limited number of echovirus types were 247 analysed. Secondly, there were too many missing back titration results to include this variable in the 248 249 mixed-effects regression analysis. Furthermore, it is possible that other unknown VNA variables, including variations in transport conditions, were not included in the analysis. Finally, as no 250 international standards exist for the included echoviruses, we had no golden standard to compare nAb 251 252 titres with. Therefore, the median nAb titre of all participating laboratories in the reference sample was considered for normalization. In conclusion, we conducted an echovirus VNA quality assessment 253 254 in different European countries. Although providing a reference serum sample can improve standardisation, comparison between sero-surveillance studies performed by different laboratories 255 256 should be done with caution and additional research is necessary to understand the effect of different VNA variables. Based on this knowledge practical recommendation to perform VNA should be further 257 258 explored.

259

260 Funding information

- 261 The study was supported by the Ministry of Health, Welfare and Sport, the Netherlands as part of the
- 262 EV surveillance programme of the National Institute for Public Health and the Environment. One of the
- 263 participating laboratories was supported by a grant from the Instituto de Salud Carlos III (PI18CIII-
- 264 00017).
- 265

#### 266 Acknowledgements

- 267 We would like to acknowledge the European Non-Polio Enterovirus Network for facilitating this study.
- 268 Furthermore, we would like to thank the Dutch National Institute for Public Health and the Environment
- 269 (RIVM), and the Departments of Medical Microbiology of Amsterdam UMC and Microvida for providing
- the materials and necessary human resources. Finally, we would like to thank Ruth Verheyen for
- 271 providing her statistical and mathematical insight. We would like to acknowledge the European Non-
- 272 Polio Enterovirus Network https://escv.eu/european-non-polio-enterovirus-network-enpen) for
- 273 facilitating this study.
- 274

#### 275 Author contributions

- K.C., K.B., K.W., J.L.M., G.K. and H.v.E.: Study design; G.K., H.v.E., E.K., I.S., M.J., L.N.G., I.G., A.S., S.D., S.B., M.C., I.T., Z.H.,
  P.S. and E.H.: Experiments; K.C. and M.K.: Data analysis; K.C., K.B., K.W., J.L.M., D.P., H.H., J.L.B. and M.B.: Conclusions
  and writing.
- 279

#### 280 Conflicts of interest

- 281 The authors declare that the research was conducted in the absence of any commercial or financial
- relationships that could be construed as a poten- tial conflict of interest.

- 284
- 285
- 286
- 287
- 288
- 289
- 290
- 291
- 292
- 293
- 294

### 295 References

- 296 1. Brouwer L, Moreni G, Wolthers KC, Pajkrt D. World-wide prevalence and genotype distribution of
- 297 Enteroviruses. Viruses 2021;13:434.
- 2. Bubba L, Broberg EK, Jasir A, Simmonds P, Harvala H, et al. Circula- tion of non-polio enteroviruses in 24 EU
- and EEA countries between 2015 and 2017: a retrospective surveillance study. Lancet Infect Dis
- 300 2020;20:350–361.
- 301 3. de Crom SCM, Rossen JWA, van Furth AM, Obihara CC. Enterovirus and parechovirus infection in children: a
   302 brief overview. Eur J Pediatr 2016;175:1023–1029.
- 303 4. Rantakallio P, Leskinen M, von Wendt L. Incidence and prognosis of central nervous system infections in a
- birth cohort of 12,000 children. Scand J Infect Dis 1986;18:287–294.

- 305 5. MeijerA, van der SandenS, Snijders BEP, Jaramillo-GutierrezG, BontL, et al. Emergence and epidemic
- 306 occurrence of enterovirus 68 respira- tory infections in The Netherlands in 2010. Virology 2012;423:49–
   307 57.
- 308 6. González-Sanz R, Casas-Alba D, Launes C, Muñoz-Almagro C, Ruiz-García MM, et al. Molecular
- 309 epidemiology of an enterovirus A71 outbreak associated with severe neurological disease, Spain, 2016.
- 310 Euro Surveill 2019;24:1800089.
- 7. Aliabadi N, Messacar K, Pastula DM, Robinson CC, Leshem E, et al. Enterovirus D68 infection in children
  with acute flaccid myelitis, Colo- rado, USA, 2014. Emerg Infect Dis 2016;22:1387–1394.
- 8. Casas-Alba D, de Sevilla MF, Valero-Rello A, Fortuny C, García-García J-J, et al. Outbreak of brainstem
- encephalitis asso- ciated with enterovirus-A71 in Catalonia, Spain (2016): a clinical observational
- 315 study in a children's reference centre in Catalonia. Clin Microbiol Infect 2017;23:874–881.
- Grapin M, Mirand A, Pinquier D, Basset A, Bendavid M, et al. Severe and fatal neonatal infections linked to
  a new variant of echovirus 11, France, July 2022 to April 2023. Euro Surveill 2023;28:22.
- a new variant of echovirus 11, France, July 2022 to April 2023. Euro Surveill 2023;28:22.
- 318 10.Fischer TK, Simmonds P, Harvala H. The importance of enterovirus surveillance in a post-polio world.
- 319 Lancet Infect Dis 2022;22:e35–e40.
- 320 11.Harvala H, Benschop KSM, Berginc N, Midgley S, Wolthers K, et al. European non-polio enterovirus
- 321 network: introduction of hospital- based surveillance network to understand the true disease burden of
- non-polio enterovirus and parechovirus infections in Europe. Micro- organisms 2021;9:1827.
- 323 12.Pons-Salort M, Grassly NC. Serotype-specific immunity explains the incidence of diseases caused by
- human enteroviruses. Science 2018;361:800–803.
- 325 13.Karelehto E, Koen G, Benschop K, van der Klis F, Pajkrt D, et al. Enterovirus D68 serosurvey: evidence
- for endemic circulation in the Netherlands, 2006 to 2016. Euro Surveill 2019;24:35.

327	14.Westerhuis B, Kolehmainen P, Benschop K, Nurminen N, Koen G, et al. Human parechovirus
328	seroprevalence in Finland and the Netherlands. J Clin Virol 2013;58:211–215.
329	15.Couderé K, van der Straten K, Brouwer L, Koen G, van Eijk H, et al. Neutralising antibodies against
330	enterovirus and parechovirus in IVIG reflect general circulation: a tool for sero-surveillance. Viruses
331	2021;13:1028.
332	16.Kamau E, Harvala H, Blomqvist S, Nguyen D, Horby P, et al. Increase in enterovirus D68 infections in
333	young children, United Kingdom, 2006-2016. Emerg Infect Dis 2019;25:1200–1203.
334	17.Zhu R, Cheng T, Yin Z, Liu D, Xu L, et al. Serological survey of neutralizing antibodies to eight major
335	enteroviruses among healthy population. Emerg Microbes Infect 2018;7:2.
336	18.Kamau E, Bessaud M, Majumdar M, Martin J, Simmonds P, et al. Estimating prevalence of Enterovirus
337	D111 in human and non- human primate populations using cross-sectional serology. J Gen Virol
338	2023;104:11.
339	19.Rosenfeld AB, Shen EQL, Melendez M, Mishra N, Lipkin WI, et al. Cross-reactive antibody responses
340	against nonpoliovirus Entero- viruses. mBio 2022;13:e0366021.
341	20.R Core Team. Vienna, Austria: R Foundation for Statistical Computing; 2020. https://www.R-
342	project.org/
343	21.Reed LJ, Muench H. A simple method of estimating fifty percent endpoints. Am J Hyg 1938;27:493–
344	497.
345	22.Nguyen D, Simmonds P, Steenhuis M, Wouters E, Desmecht D, et al. SARS-CoV-2 neutralising antibody
346	testing in Europe: towards harmonisation of neutralising antibody titres for better use of
347	convalescent plasma and comparability of trial data. Euro Surveill 2021;26:27.
348	23.Khoury DS, Wheatley AK, Ramuta MD, Reynaldi A, Cromer D, et al. Measuring immunity to SARS-

349	CoV-2 infection: comparing assays a	nd animal models. Nat Rev Immunol 2020;20:727–738
-----	-------------------------------------	---

- 350 24.Manenti A, Molesti E, Maggetti M, Torelli A, Lapini G, et al. The theory and practice of the viral dose
- 351 in neutralization assay: insights on SARS-CoV-2 "doublethink" effect. J Virol Methods
- 352 2021;297:114261.
- 353 25.Cooper G. Report of the WHO collaborative study to establish the 1st WHO International Standard for
- 354 anti-EV71 serum (Human). WHO/BS/2015.2267. 2018.
- 26.Crawt L. Report on the WHO collaborative study to establish the 1st International Standard for Sabin
- 356 inactivated polio vaccine (sIPV). WHO/BS/2018.2338; 2018

Table 1. VNA characteristics of different participating laboratories. Origin of cell lines:

RD=Rhabdomyosarcoma, HT29=Human B cell lymphoma cell line, A549=Alveolar carcinoma epithelial cell line, Caco-2=Colorectal adenocarcinoma cell line, LLC- MK=Rhesus Monkey Kidney Epithelial Cell line 2. Abbreviations: TCID<sub>50</sub>=50% tissue culture infectious dose,

CPE=Cytopathogenic effect, nAb=Neutralising antibody

Participa ting laborato ry	Cell line	Mono- layer (M)/ Suspens ion(S)	No. of cells/well	TCID <sub>50</sub> virus/ well adde d	Serum sampl e inacti vation	Virus /sampl e co- incuba tion before adding cells (hours)	Assay read-out	Method of nAb titre calculati on	Incu batio n lengt h (days)	Chloroform used
1	HT29, RD, A549	S	Not counted	100	Yes	1	CPE	Reed Muench	7	Yes (E18)
2	RD	М	18000	100	No	3	CPE	Reed Muench	6	No
3	RD	Μ	15000	100	No	2	CPE	Last dilution with 100% neutralisati on	3	No
4	RD	Σ	20000	100	No	1	CPE	Improve d Kar ber Met hod	5	No
5	RD	S	15000	100	Yes	1	CPE	Reed Muench	3	No
6	RD	М	25000	100	No	3	CPE	Reed Muench	7	No
7	RD	М	25000	200	No	3	Crystal violet	Reed Muench	3	No
8	RD	S	15000	100	No	3	Crystal violet	Reed Muench	6	No
9	Caco-2	S	Not counted	100	Yes	3	CPE		5	No
10	RD	S	14400	100	Yes	1.5	CPE	Last dilution with 100% neutralisatio n	6	No
11	RD	М	Not counted	100	No	3	CPE	Last dilution with 100% neutralisatio n	4	No
12	RD, LLC- MK2	S	15000 (10000 LLC- MK2)	100	Yes	1	CPE	Last dilution with 100% neutralisatio n	RD 3, LLC- MK2 4	No

Table 2. Geometric coefficients of variation (GCV) of nAb titres determined by ten laboratories in IVIG dilutions 1:1 (S1), 1:4 (S2), 1:16 (S3) before and after calibration

		E30 G2	E30 G6			
	Measured GCV (%)	Calibrate GCV (%)	Measured GCV (%)	Calibrated GCV (%)		
S1	10.3	9.7	15.6	15.8		
S2	24.5	11.6	50.6	16.6		
S3	60	17	37.8	31.5		

Fig. 1. Neutralising Ab titres by participating laboratories according to their proper SOP in (a) S1 against echovirus types E1, E18, E30 G2, E30 G6 and E6 (b) dilution series S1–S4 against E6.



Fig. 2. Comparison between measured and calibrated nAb titres against E30 G2 (a) and E30 G6



